NAVODILO ZA UPORABO

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DIRECTION FOR USE

1. Usage

Thermo knife REF 0605 is electrical device for waxing up in dental laboratory.

2. Advantages

Simple use, efficient, always ready to be used.

3. Technical data

Connection	230V,50/60Hz
Fuse	0,16A
Power taken away at max. temperature	8W
Heater	voltage 6V/5W, cable 1,2m
Temperature range:	30 - 250°C not level regulation
Function signalization:	Lightening switch
Outer measures:	w=120mm, d=80 mm, h=60 mm
Weight:	0,8 kg
Cable with protective plug	1,5 m

4. Description

REF 0605 Thermo knife is compact, efficient device for waxing up by hand with hot point. Level regulated temperature is not possible and this provide conditions for waxing up tension free. The range of temperature regulation is 30°C. Cork that isolate the holder make using the thermo knife safe.

4.1. Before first use

When you take new device from the box make sure it is not damaged from the transport. If you notice any damage return the device immediately. Put the device on desired place and connected it with the plug on heating container. Make sure that the net voltage is in accordance with data on the board of device. If it is appropriate than connect the device on the net.

4.2. Instruction

4.2.1. Direction

It is very easy to change wax-up instruments. You can do it by hand when the thermo knife is cold and turned off. During work you can change the instruments by special tongs resistant to temperature 250°C. **Be careful on burns!** Turn on the device with switch. It is ready when there is a green light on the switch. The desired working temperature is placed with swivel knob. Heating point is than prepared for use.

4.2.2. Mistakes removing

If wax up instruments do not heat , if electronic do not work and if even when changing the heating pin the device does not give desired results, please contact authorized service that will find and remove the mistake.

4.2.3. Maintenance

Maintain clean device. Repairing of device must always perform authorized service. Interdent is not responsible for mistakes made by unauthorized services. Guarantee for device is one year and does not valid in case of using the device in other purposes that are described in this directions for use and in case of intervention in device by unauthorized service.

5. Safety precaution

Burn risk at direct contact of hot point for wax up or hot zone of heating pin. When heating pin is not used, it must be lay down on holder to prevent overheating and possible fire. The intervention in device can be made by qualified person. Do not use device on people.

6. Content of supply

REF 0605 Thermo knife	1 pcs
REF 0617 Cork holder	1 pcs
REF 0606 Wax up instrument knife	1 pcs
REF 0607 Wax up instrument awl	1 pcs
REF 0608 Wax up instrument point thin 0,5 mm	1 pcs
REF Wax up instrument point middle 1,0 mm	1 pcs

7. Accessory

REF 0618 Cork holder
REF 0610 Wax up instrument knife, small
REF 0611 Wax up instrument knife, small, transverse
REF 0612 Wax up instrument, with hollow
REF 0613 Wax up instrument, opener, left, 0,5 mm
REF 0614 Wax up instrument, opener, left, 1,0 mm
REF 0615 Wax up instrument point, extra thin- hart gold plated 0,3 mm

8. Warranty

Our recommendations and instructions are based on our own experiences. Using of device in any other purposes as described above is made on own responsibility. Quality of device is tested and assured.

SI - 3000 CELJE - Opekarniška 26 - Tel.: +386 (0)3 425-62-00 - Fax: +386 (0)3 425-62-10
E-mail: info@interdent.cc - http://www.interdent.cc

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